Read this instruction sheet thoroughly before use!

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412-73791
0318D6

Intended use

Ketone Body Calibrator • 300 is designed to be used with Autokit Total Ketone Bodies and Autokit 3-HB for the determination of total ketone bodies (T-KB) and 3-hydroxybutyrate (3-HB) in serum and plasma.

Reagents and ingredients

Ketone Body CAL: 4 x 5 mL
Contains: 300 µmol/L 3-HB and 0.025 % sodium azide

Determination of assigned values

3-HB is specifically oxidized by 3-hydroxybutyrate dehydrogenase in the presence of -nicotinamide adenine dinucleotide, oxidized form (NAD) to produce acetoacetate and -nicotinamide adenine dinucleotide, reduced form (NADH). The 3-HB concentration is determined by measuring the amount of produced NADH, which is calculated with the molar extinction coefficient 340 nm.

Preparation of reagent

Use the calibrator as supplied.

After opening the bottle, it is recommended to use the calibrator immediately. If the calibrator is stored after opening, tightly cap the bottle and store it at 2–10 °C.

Warnings and Precautions

- For *in-vitro* diagnostic use.
- For professional use.
- Do not use the reagents described above in any procedures other than those described herein. Read instruction sheets of Autokit Total Ketone Bodies and Autokit 3-HB.
- Do not use the containers for any purpose other than those described herein.
- Do not use the reagents past the expiry date.
- After opening the reagents, it is recommended to use them immediately. When the opened reagents are stored, cap the bottles and keep them under the specified conditions.
- If the reagents come in contact with mouth, eye or skin, wash off immediately with plenty of water. Consult a physician if necessary.
- When discarding the reagents, dispose of them according to local or national regulations.
- This Calibrator Set contains sodium azide as a stabilizer. Sodium azide may react with lead or copper plumbing to form explosive compounds. Even though the reagents contain minute quantities of sodium azide, drains should be flushed well with a large amount of water, when discarding the reagents.

Ordering information

Code No.	Product	Package
412-73791	Ketone Body Calibrator • 300 (3-Hydroxybutyrate: 300 µmol/L)	CAL: 4 x 5 mL

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