



CRP-HS II LT



METHOD: Turbidimetric Immunoassay

Ordering information

Reagent	Package	Code No.
CRP-HS II LT R1 (Buffer)	R1: 2 x 40 mL	462-39493
CRP-HS II LT R2 (Latex Reagent)	R2: 2 x 40 mL	468-39593
CRP-HS LT Calibrator Set HO	CAL: 5 conc. x 2 mL	460-25801

Measuring Temperature 37°C

AU2700 / AU5400 (BECKMAN COULTER)

Parameter Setting	Multi-Point Calibration
Reagent	CRP-HS II LT
Specific test parameters	
Test No. / Name	* / CRP
Sample	Serum
Operation	YES
Sample Vol / Dil	2.5 / 0
Reagent 1 Vol / Dil	150 / 0
Reagent 2 Vol / Dil	150 / 0
Wave (Main / Sub)	600 /
Dilution	1
Method	RATE
Reaction	+
Point 1 (1st / Lst)	13 / 27
Point 2 (1st / Lst)	
Linearity (1st / 2nd)	100%
No-Lag-Time	NO
Min OD / Max OD (L / H)	0.0000 / 2.5000
Reagent OD (1stL / 1stH)	
Limit (LstL / LstH)	
Dynamic range (L / H)	
Correction Factor (A / B)	1.000000 / 0.000000
Calibration Parameters	
Calibration Type	6AB
Formula	Spline
Cuvettes / Counts	/ 4
Process	CONC
Calibrator	Cal # - OD - Conc- Factor
Point 1	* - [] - *1 -
Point 2	* - [] - *1 -
Point 3	* - [] - *1 -
Point 4	* - [] - *1 -
Point 5	* - [] - *1 -
Point 6	* - [] - *1 -
Point 7	
1 point cal point	
MB type Factor	

* : To be inserted by the user.

0311D1/vk

*1: Assigned value of the calibrator (mg/L or mg/dL).